ENA  **post**

Perno endocanalare in fibra di vetro
Glass Fiber Post
Glasfiber-Wurzelstift

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Glass fiber post - ENAPOST

Perfect aesthetic and functional integration

- Fibers are inserted in the same resin matrix of the microhybrid composite for restoration
- Aesthetic: dentine colour and fluorescent
- Functional: same elasticity of dentine
- Perfect cohesion fibers-composite
- Radiopaque

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
The coloured rings can be used to measure the depth of the canal, positioning it where the post must be cut. Note: the rings must be removed before cementation.

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Extracted tooth endodontically treated. UV light shows again a good fluorescence.

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Tooth after 35% hydrogen peroxide bleaching treatment (EnaWhite). UV light shows lower fluorescence.

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Composite restoration with EnaPost Fiber Post. UV light shows a better fluorescence than before restoration.

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Adhesive relation: dentine - enamel
dentine adhesive - composite
cement - post.

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Section of restored tooth using Enapost conicity 10% and Enamel plus HFO (picture under UV light in the middle, transilluminated section on the right).

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Section of restored tooth using Enapost conicity 2% and Enamel plus HFO (transilluminated section on the right).

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Restoration of endodontically treated anterior tooth with 2% conicity EnaPost and Enamel Plus HFO

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
Fig. 50.
 Rx finale.
Final X-Ray.
Abschluss-Röntgenbild.

(GB) Instructions for Use

A. Complete the endodontic therapy, and prepare the tooth for the restoration by evaluating if a post is necessary, based on residual quantity of intact tooth. The loss of one or both interproximal walls in anterior teeth, of one or both mesial and distal cuspid in posteriors, suggest you to use a post, in order to better stabilize the restoration.

B. Make an X-ray to determine the appropriate diameter and depth of post space preparation. The post diameter should be slightly lower than the endodontic preparation. The post should be placed at same radicular depth as the height of dentinal body (for composite aesthetic restoration) or of post and core. **Use of rubber dam is recommend.**

C. Remove gutta-percha to preplanned depth with a Gates-Glidden drill or Peeso reamer. Radiographic verification is recommended.

D. Select the EnaPost Drill diameter which corresponds to the last Gates-Glidden drill to begin preparing the post space.

Use Drills with a slow-speed contra-angle with water spray in order to avoid tissues overheating. The drill must be kept in continuous clockwise rotation until it has been completely removed from the tooth. This will minimize the risk of the drill seizing in the post space. Remove any tooth debris from canal by irrigating post space with water spray. Sequentially step up to the next larger **EnaPost** Drill until the preplanned diameter and depth are achieved.

E. Select the **EnaPost** that corresponds to the last drill used to prepare the post space and insert post into post space. For a better fitting shorten the post as necessary from either the apical or occlusal end (as clinical judgement dictates) with a Diamond separating disc, wetting the post with water to avoid fibers overheating.

F. Remove the post and gently clean it with an alcohol wipe and apply on the surface a mixture of **EnaBond** with **Enabond Catalyst** (one drop+one drop).

G. Etch the cavity (coronal and radicular preparation) with EnaEtch 37% phosphoric acid for 2 minutes. For ideal adhesion, before etching, a microblasting on cavity surfaces is recommended in order to clean and eliminate endodontic material debris.

H. Wash accurately the canal with a syringe to remove completely the acid. Suck water and dry the canal with paper point; in order to maintain wet the dentine and avoid collagen collapse do not dry with air.

I. Apply in the cavity and in the canal the mixture of **EnaBond** and **EnaBond Catalyst**.

Ricostruzione di perno moncone con EnaPost (conicità 10%)
Restoration of Post and core with EnaPost (conicity 10%)
Restauration durch Stift-/Stumpf Aufbau mit Enapost (10% Konizität)

Fig. 51.
 Moncone protesico gravemente compromesso
Severly compromised prothetic core
 Äusserst geschwächter prothetischer Stumpf

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The adhesive should be wiped on the surface with a disposable microbrush or with a paper point; dry with air to eliminate water and solvent residuals. Insert the post to check the canal and better push adhesive in dentine tubules.

J. Apply a disposable automixing tip and extract **EnaCem** dual composite cement from the syringe: the two components will be automatically mixed. Place the cement directly into the canal filling it from the bottom of the cavity to the surface using a disposable endodontic applicator (Hawe-Neos); do not use any lentulo or rotating instrument. Apply some cement on the post surface and insert slowly the post to full depth, allowing excess cement to vent. Find the right position and, maintaining the post in position, wipe away any excess cement.

Note: Enacem working time 4 min. Complete setting time 7-9 min.

K. Light cure for 60 seconds and proceed to the restoration using **Enamel Plus HFO** Micro-Hybrid composite, both for aesthetic direct restoration and post and core, that could be immediately prepared once the restorative procedure has been completed.

Cleaning and Sterilization

Drills: Prior to each use, remove debris ultrasonically or by hand scrubbing with a brush; thoroughly rinse and dry. Place drills, into an autoclave bag and seal. Sterilize at 132 C in compliance to a system validated by the dentist. Sterilizing bags should not touch the walls of the autoclave.

Posts: Posts may be sterilized **only** by cold sterilant

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MICERIUM S.p.A. Via G. Marconi, 83 -16036 Avegno (GE) Italy